-ebro-
a xylem brand

**pH-Meter
pH-meter
pH-mètre**

PHT 810

Inhalt

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
Table of Contents

Overview	5
What is the meaning of the term.....	7
Safety instructions	9
Unpacking / scope of delivery	10
Operation	15
Switching the device on/off.....	15
Performing measurements	17
Holding the current measurement	19
Changing the device settings	21
Changing the sensor	24
Calibrating the device	27
Preparing for calibration	27
Calibration.....	29
Changing the battery	35
What to do, if...	38
Maintenance and disposal	41
Cleaning.....	41
Disposal	41
Technical data	43
Approvals	45
Operation	45

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Overview

The PHT 810 can be used to measure the pH value of a liquid.

Various sensors are available for different applications.

Device elements:

- 1 Display (LCD)
- 2 ON/HOLD key (on / off / store displayed measured value)
- 3 Key ▲
- 4 CAL key
- 5 Key ▼
- 6 BNC connector for probe

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What is the meaning of the term...

...BNC connector

A BNC connector is a coaxial connector with a so-called bayonet nut. Coaxial means that an inner conductor is surrounded by a hose-shaped outer conductor, which has a shielding effect due to its shape. A bayonet nut operates according to the principle “put on - turn slightly (usually approx. 90°) – connected”.

...pH-value

The pH value is a measurement for the concentration of hydrogen ions in a solution. The value is stated as a negative base 10 logarithm ($\text{pH} = -\log_{10} a(\text{H}^+)$). Acidic solutions have pH values < 7 ; alkaline solutions have pH values > 7 . pH comes from “pondus Hydrogenii” (lat.: pondus = weight, Hydrogenium = hydrogen).

...Buffer solution

A buffer solution absorbs contaminants with acids or bases without a material change in the pH value. Therefore the pH value of a buffer solution ideally remains constant, even when an acid or base is added (e.g. through immersion of a sensor coated with an acid or base).

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Safety instructions

Do not use the device in explosive areas.

Danger of death!

Do not use the device in an environment hotter than 60 °C.

The lithium battery may explode!

Please observe the following rules in order to avoid either damaging the device or making measurement errors:

- Protect the device housing from direct contact with water.
- Always use the measuring probe according to the instructions.

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Unpacking / scope of delivery

Check the device packaging and contents to ensure it is complete and undamaged. You must also check that the contents of the packaging match your order.

Scope of delivery

The delivery includes these parts:

- Device PHT 810 with sensor
- This user manual

The delivery may also include various accessories. If you have any reason for complaint, do not hesitate to contact us. You will find our contact data on the back of this manual.

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Possible accessory items

Description	Identification
Case for PHT 810	AT100
Electrode for measurements in water, galvanic...	AT200
Electrode for measurements in chemistry, pharmaceuticals	AT201
Electrode for cut-in measurements in foods	AT202
Electrode for cut-in measurements in foods, with plastic shaft	AT204
Electrode for cut-in measurements in foods, with plastic shaft and grip	AT206
Punch	AT300
Buffer solution pH4 (50ml)	AT400
Buffer solution pH7 (50ml)	AT401
KCL solution (50ml)	AT405
Electrode cleaner	AT410
Drop protection	AG140

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Operation

Switching the device on/off

- To turn the device on, press the ON/HOLD key.

After a few seconds (self-test) the device is ready to take measurements and already displays the current pH value.

A small battery symbol that shows the remaining charge in the device battery (see page 35) appears at the top left of the display.

- To turn the device off, press the ON/HOLD key until “OFF” appears on the display (several seconds).

If you forget to turn off the device it will turn off automatically – provided this setting has not been modified (see page 23) – after two hours of operation to conserve the battery.

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
Performing measurements

In order to achieve the highest level of measurement accuracy, you need to determine the temperature of the substance to be measured and set this on the device each time before you measure the pH value. To perform measurement in solid or semi-solid media and to avoid damage of the PH probe, use the accessory drill AT300.

Setting the temperature

- Briefly press the CAL key.
- Press the keys or until the actual temperature of the substance to be measured is shown on the display.
- Briefly press the CAL key.
- Press the ON/HOLD key.

Carrying out the measuring process

- Insert the sensor into the substance to be measured and read the displayed pH value.

If the display shows the message “HI” or “LO”, please see page 39, “What to do, if...”.

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
Holding the current measurement

You can hold the current reading on the display in order to better recognise changes in the measured value.

- In order to hold the current reading, briefly press the ON/HOLD key one times.

A display line will appear below the current measured value. This line shows the indicator "HOLD" and the held reading to the right of it.

The measured value in the top display line is continuously updated (measuring cycle rate see page 23, RATE).

The held reading is displayed until you press the ON/HOLD key again.

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Changing the device settings

You can set:

- Whether the device displays the temperature in °Celsius or °Fahrenheit,
- Whether the automatic device switching should be activated or not,
- How often the measured value is to be updated (every 1 to 15 seconds).

You must call up the User-Menu to change the device setting. For this, proceed as follows:

- Make sure that the device is switched off
- Press and hold the ON/HOLD button until the indicator "PHT" appears in the display
- Press and hold the CAL key until the indicator "USER" appears

The device is now in configuration mode and displays the User Menu. This consists of the menus UNIT, SHDWN and RATE.

- Page to the next menu point using the ▼ key. Page back to the previous menu point using the ▲ key.
- Change the value with the CAL key. Changed values will be saved immediately.

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
UNIT

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SHDWN

Automatic shutdown after 2 hours of operation

ON = shutdown activated
OFF = shutdown de-activated

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RATE

Measuring cycle rate (1...15 seconds); above 15 the scale reverts to 1

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Modified values are automatically stored right away.

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- To leave the User-Menu, switch off the device by pressing the ON//HOLD button.

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Changing the sensor

The sensor is connected to the device with a BNC connector.

To change the sensor, proceed as follows:

- Turn the sensor on the BNC connector by 90° counter-clockwise while slightly pressing on the sensor in the direction of the device.
- Pull the sensor on the connector out downwards.
Do not pull on the sensor!
- Insert a new sensor.
- Turn the sensor and/or the extension cable on the BNC connector 90° clockwise.
The BNC connector must lock in place!
- Ensure the sensor connection is firmly seated.

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Calibrating the device

In order to maintain a high level of measurement accuracy, the device must be calibrated every time the sensor is changed and after every three months of use at the latest. During calibration the gradient of the device is compared to one, two, or three buffer solutions that have a defined pH value.

The larger the number of buffer solutions used and therefore the more points on the gradient are calibrated, the higher the subsequent measurement accuracy of the device will be.

The device has a special calibration mode that guides you through the calibration.

Preparing for calibration

To calibrate the device, you require at least one and at the most three buffer solutions with different pH values, which must be exactly 4.01, 7.00, or 10.01. The temperatures of these buffer solutions must be close to the same.

- If required, measure the temperature of the buffer solutions and ensure that the temperatures are close to the same during use of several buffer solutions.

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
NOTE:

Material temperature differences between the buffer solutions reduce the quality of the calibration and therefore the subsequent measurement accuracy of the device.

- Remove the protective cap from the sensor and rinse the sensor for approx. 10 minutes in a pH neutral liquid or in tap water.

Calibration

- To carry out the calibration, briefly press the CAL key. This starts the calibration mode of the device.

Calibration can be discontinued at any time. To do so, briefly press the ON/HOLD key.

Discontinuing calibration or an “ERROR” while determining the calibration values (see below) leads to lower measurement accuracy.

If the device does not receive any key inputs for several minutes while in calibration mode, it automatically switches back to measurement mode.

- First, enter the temperature of the buffer solution. To do so, press the ▼ or ▲ keys until the actual temperature of the buffer solution(s) is shown on the display.
- Press the CAL key.

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

- Next, specify if you want to calibrate the device gradient using one, two, or three buffer solutions. To do so, press the ▼ or ▲ , key until “1PC”, “2PC”, or “3PC” is shown.
- Press the CAL key. The display shows “CAL 1”.
- Insert the sensor into the first buffer solution and then press the CAL key again.

The device now determines the first calibration value. If the calibration value is valid, “OK” appears on the display after some time; otherwise “ERROR” appears (see note on page 27). In the latter case, repeat the process with another buffer solution.

If you previously selected “2PC” or “3PC”, “CAL 2” is now shown on the display. Otherwise, calibration is complete and the device switches back to measurement mode.

- If “CAL 2” is shown on the display, insert the sensor into the second buffer solution and press the CAL key.

The order of the buffer solutions with the various pH values is random.

The device now determines the second calibration value. If the calibration value is valid, “OK” appears on the display after some time; otherwise “ERROR” appears (see note on page 27). In the latter case, the device requests the second buffer solution (“CAL 2”) again.

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
If applicable, use a buffer solution with a pH value that has not yet been used during this calibration session.

If you previously selected “3PC”, “CAL 3” is now shown on the display. Otherwise, calibration is complete and the device switches back to measurement mode.

- If “CAL 3” is shown on the display, insert the sensor into the third buffer solution and press the CAL key.

The device now determines the third calibration value. If the calibration value is valid, “OK” appears on the display after some time; otherwise “ERROR” appears (see note on page 27). In the latter case, the device requests the third buffer solution (“CAL 3”) again.

If applicable, use a buffer solution with a pH value that has not yet been used during this calibration session.

If the device has determined a valid third calibration value, calibration is complete and the device switches back to measurement mode.

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


Changing the battery

The battery symbol in the display shows the battery condition.


To change the lithium battery, the device must be opened. Proceed as follows:

Use suitable precautions such as the use of a grounded wrist strap in order to prevent electrostatic discharge while changing the battery! Electrostatic discharges can destroy the device!

- First, carefully remove the two rubber stoppers (1) on the lower side of the device, e.g. using forceps.
- Now completely remove the two visible Phillips-head screws (2) with a suitable screwdriver (PZ1).
- Hold the bottom of the device (3) with one hand and carefully pull the upper part (4) off upwards with the other hand until the circuit board is exposed.

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
Now you can see the battery mount (7) with the battery (6) on the circuit board (5).

- Pull the spent battery from the battery mount in the direction of the arrow.
- Take the new battery with fat-free fingers and insert it into the battery mount. The plus sign on the battery must face upwards so it is visible.
- Reassemble the upper and lower parts of the housing.
- Tighten the two screws (only with moderate force, 0.4 Nm) and replace the rubber stoppers (the slightly beveled edges facing out).

Dispose of the spent battery in an environmentally sound manner!

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Possible cause	Remedy
Measurement range exceeded	Observe Measurement range
Sensor defective	Call service
Below measurement range	Observe measurement range
Short-circuit of the sensor	Call service

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Maintenance and disposal

Cleaning

Clean the device with a slightly dampened cloth.

Never use solvents (such as acetone) for cleaning because these may attack the plastic.

Disposal

If the device becomes no longer fit for purpose, it must be disposed of in a suitable, environmentally-friendly manner.

Never dispose of the device in the domestic garbage. Instead please return it to us. We will take care of disposal in an environmentally sound manner.

Dispose of the battery at a designated collection point.

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Technical data

Measurement range	pH 0...14
Resolution of the measurement signal	pH 0,01
Measurement accuracy of the electronics	pH $\pm 0,03$
Temperature of the substance being measured	
Adjustment range	0...120 °C (32...248 °F)
Increment	0,5 °C (0,5 °F)
Operating temperature	0...+50 °C (32... 122 °F)
Category temperature	-10...+60 °C (14... 140 °F)
Storage temperature	-25...+60 °C (-13 ... 140 °F)
Refresh rate	1...15 seconds, adjustable
Automatic shutdown	automatically after 2 hours, can be deactivated
Protection class	IP 67
Dimensions	115 x 54 x 22 mm (L x W x H)
Housing material	ABS
Weight approx.	90 g
Sensor connector	BNC
Battery	Lithium 3 V / 1 Ah
Battery life	up to 5 years

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Approvals

According to the certificate of conformity, this product complies with the EMC directive.

Information, operating instructions and declarations of conformity can be found at www.ebro.com

Operation

Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) This device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications of the device could void the user's authority to operate the equipment.

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Service-Adresse / Service Address / Adress du Service

Xylem Analytics Germany Sales GmbH & Co. KG
ebro
Am Achalaich 11
82362 Weilheim
Germany

Phone: +49.(0)841.954.78.0
Fax: +49.(0)841.954.78.80
Internet: www.ebro.com
E-Mail: ebro@xylem.com

.....

Hersteller / Producer / Fabricant

Xylem Analytics Germany GmbH
Am Achalaich 11
82362 Weilheim
Germany