



Non-contact

±0.2°C

Accuracy

250ms

Response time



Data hold

1cm~5cm

Test distance



High temperature alarm



UT30H



Specifications	UT30H	UT308H
<b>Certificates</b>	CE <sub>0413</sub> , FDA, UKCA	
Measurement range (object)	0°C~100°C	0°C~100°C
Accuracy (object)	±1°C	±1°C
Measurement range (body)	32°C~45°C	32°C~45°C
Accuracy (body)	±0.2°C	±0.2°C
Resolution	0.1°C	0.1°C
Data storage		32 groups
Optimum measuring distance	1~5cm	
Unit conversion (°C/°F)	√	
Measurement alarm	Audible alarm for >37.2	
Power on self test	√	
Auto shutdown	√	
Backlight	√	
Operating environment	15°C~35°C, <85%RH	
Transport and storage environment	-20°C~60°C, <93%RH	
Operating atmosphere pressure	700hPa-1060hPa	

Characteristics	
Standard accessories	Batteries
Power	UT30H: 1.5V AAA x 2; UT308H: 1.5V AA x 2
Product size	UT30H: 141 x 90 x 36mm; UT308H: 157 x 99 x 48mm
Product net weight	UT30H: 98.7g; UT308H: 140.8g
Standard individual packing	Gift box, blister, English manual
Standard quantity per carton	20pcs
Standard carton measurement	UT30H: 310 x 270 x 210mm; UT308H: 405 x 320 x 270mm
Standard carton gross weight	UT30H: 4.5kg; UT308H: 6.6kg